

**WATER TREATMENT SOLUTION****DESCRIPTION**

**DYNASPHER WA907-WT** IS A MACRO POROUS WEAK ANION EXCHANGE RESIN HAVING TERTIARY AMINE GROUP. IT HAS EXCELLENT MECHANIC CHARACTERISTICS, THANKS ITS POROUS STRUCTURE ADSORB AND DESORBS ORGANIC SUBSTANCE WITH STRONG FACILITY.

**DYNASPHER WA907-WT** SHOWS HIGH REGENERATIVE EFFICIENCY, IT IS SUITABLE IN CO-CURRENT OR COUNTER-CURRENT DEMINERALIZATION SYSTEM, FLOATING BED, LAYERED BED.

**SYSTEM DESIGN**

Co - current / Counter current / Floating bed / Blocked bed

**PRINCIPAL APPLICATIONS**

- Fruit juices demineralization
- Fruit juices deacidification
- Liquid sugar
- Milk whey
- Pharmaceutical
- Nutraceutical
- Metallurgical
- water

**REGULATORY**

- F.D.A. – CFR 21 – 173.25
- Codes Alimentarius – Inventory of Processing Aids – CAC/MISC3
- European Resolution AP (97) – 1 regarding the TOC (Total Organic Carbon) released according AFNOR method (method T90 – 601)

**TYPICAL PACKAGING**

- 1 ft<sup>3</sup> Sack
- 25 lt Sack
- 5 ft<sup>3</sup> Drum (Fiber)
- 1 m<sup>3</sup> Supersack
- 42 ft<sup>3</sup> Supersack



PRODUCT INFORMATION  
**DYNASPHER WA907-WT**  
MACROPOROUS WEAK ANION RESIN

## TYPICAL CHARACTERISTICS

### PHYSICAL CHARACTERISTIC

Copolymer	Styrene-divinylbenzene
Matrix	Macroporous
Type	Weak anion
Functional Group	Tertiary/Quaternary amine
Physical Form	Beige opaque spherical beads

### CHEMICAL CHARACTERISTICS

Ionic Form as Shipped	FB
Total Exchange Capacity	≥ 1.6 eq/lit
Water Retention	50.0 - 58.0 %

### PARTICLE SIZE

Particle size range	0.425 mm - 1,2 mm
Uniformity Coefficient	≤ 1.5
< 300 µm	≤ 1.0 %
> 1180 µm	≤ 1.0 %

### STABILITY

Whole Uncracked Beads	≥ 98 %
Swelling	Cl <sup>-</sup> → OH <sup>-</sup> - 25% max

### DENSITY

Particle Density	1.060 - 1.100 g / ml
Shipping Weight	650 - 750 g / lit

For additional size in formation, please refer to the our Technical Dept.

### SUGGESTED OPERATING CONDICTIONS

Termal stability	70 °C FB (158 °F) – 100 °C (212 °F)
Chemical stability	0 - 9
Service Cycle	0 - 14
Minimum bed depth	800 mm

For additional praticle size information regarding recommended minimum bed depth, operating conditions, and regeneration conditions for Layered or Mixed bed, please refer to our tecnical dept.

# PRODUCT INFORMATION

## DYNASPHER WA907-WT

### MACROPOROUS WEAK ANION RESIN

## HYDRAULIC CHARACTERISTICS

Estimated bed expansion of DYNASPHER WA907-WT as a function of backwash flowrate and temperature is shown in figure 1.

Estimated pressure drop for DYNASPHER WA907-WT as a function of service flowrate and temperature is shown in figure 2.

These pressure drop expectations are valid at the start of the service run with clean water and well – classified bed.

Figure 1: Back wash expansion  
Temperature = 5 °C - 40 °C (41 °F - 104 °F)





Figure 2: Pressure Drop  
Temperature = 5 °C - 40 °C (41 °F - 104 °F)



## CUSTOMER NOTICE

### STORAGE

It is recommended to store ion exchange resins at temperatures above the freezing point of water under roof in dry conditions without exposure to direct sunlight. If resin should become frozen, it should not be mechanically handled and left to thaw out gradually at ambient temperature. It must be completely thawed before handling or use. No attempt should be made to accelerate the thawing process.

### DISPOSAL

In the European Community ion exchange resins have to be disposed, according to the European waste nomenclature which can be accessed on the internet – site of the European Union.

### TOXICITY

The safety data sheet must be observed. It contains additional data on product description, transport, storage, handling, safety and ecology.

### WARNING

Oxidizing agents such as nitric acid attack organic ion exchange resins under certain conditions. This could lead to anything from slight resin degradation to a violent exothermic reaction (explosion). Before using strong oxidizing agents, consult sources knowledgeable in handling such materials.